

## PROTEIN G COATED SURFACES

The Biomat product is a 96 well coated microplate with recombinant Protein G and a protein to block non-specific binding sites and to maintain stable activity.

Protein G specifically binds the Fc region of immunoglobulins of many mammalian species ( see table 1 ), with an orientation that allows the F(ab)<sub>2</sub> binding sites to be freely available for efficient binding to epitope. When coated onto microplates, the Protein G can securely capture IgG applied directly or as antigen/antibody complexes.

### *Example of applications:*

- *specific and sterically oriented bond of IgG*
- *separation of IgG from other immunoglobulins*
- *separation of antigen-antibodies complexes*
- *isolation and analysis of fusion proteins*

### **Product specifications**

#### **Components**

Individually pouched 96-well microplates, configured in 12 removable 8-well strips

#### **Coating**

Recombinant Protein G ( mol. weight 26.1 kDa ), from Streptococcus sp., expressed in E.coli, is coated using 100 µl/well. The strips are post-coated ( blocked ) for low non specific binding and long-term stability.

#### **Binding capacity**

Microplate was saturated with human IgG at a concentration of 8.0 µg/ml ( 800 ng/well ) in an ELISA format using Streptavidin-HRP as detector and TMB as substrate ( see Figure 1 for data and experiment details ).

The Biomat Protein G microplate shows a nominal **binding capacity of ~ 5.3 pmol IgG/well**

#### **Sensitivity**

Biotinylated human IgG was detected at a concentration significantly above background in an ELISA format using streptavidin-HRP as detector and TMB as substrate ( see Figure 2 for data and experiment details ).

The Biomat Protein G microplate shows a **sensitivity of 0.056 ng/well of human IgG**.

#### **Uniformity**

Microplates show a **CV% less than 5** when used as a catcher of biotinylated human IgG in an ELISA format using streptavidin-HRP as detector and TMB as substrate.

#### **Storage and Stability**

The microplates, if unopened, are stable refrigerated until the expiration date printed on the label. If opened, store in closed pouch with dessicant and use within the expiration date.

# *biomat*

Table 1. Binding affinities of recombinant Protein A and G for Immunoglobulin binding domains

( The table below gives an overview of binding strengths of protein A and G to different species and subclasses. S : strong binding; M : medium binding; W: weak binding; N : no binding )

| <b>Species</b> | <b>Ig Subclass</b> | <b>Protein A</b> | <b>Protein G</b> |
|----------------|--------------------|------------------|------------------|
| Human          | Total Ig           | S                | S                |
|                | IgG1,IgG2,IgG4     | S                | S                |
|                | IgG3               | W                | S                |
|                | IgD                | W                | N                |
|                | IgA                | W                | N                |
|                | IgE                | W                | N                |
|                | IgM                | W                | N                |
| Mouse          | Total Ig           | S                | S                |
|                | IgG1               | W                | M                |
|                | IgG2a, IgG2b,IgG3  | S                | S                |
|                | IgM                | N                | N                |
| Rabbit         | IgG                | S                | S                |
| Rat            | IgG                | N                | W-S              |
| Goat           | IgG                | W-M              | M-S              |
| Sheep          | IgG                | W-M              | M-S              |
| Chicken        | IgG                | N                | W                |
| Guinea Pig     | IgG                | S                | W-M              |
| Hamster        | IgG                | W                | M                |
| Horse          | IgG                | W                | S                |
| Pig            | IgG                | S                | W-M              |
| Bovine         | IgG                | M                | S                |
| Dog            | IgG                | S                | W-M              |
| Cat            | IgG                | S                | W                |

## TECHNICAL NOTES N. 26 Binding capacity test


1. Add 100  $\mu$ l of different concentrations of human IgG ( from 0.5 to 10  $\mu$ g/ml ) mixed with a constant amount of human IgG biotinylated ( 0.01  $\mu$ g/ml ) to the wells of Protein G coated plate and incubate for 30 minutes at room temperature.
2. Empty the wells and wash with 0.1 M PBS pH 7.2,0.05% Tween<sup>®</sup> 20 four times
3. Add 100  $\mu$ l /well of Streptavidin-HRP ( BioSpa product code SB01-61, diluted 1:20.000 ) and incubate for 30 minutes at room temperature.
4. Empty the wells and wash with 0.1 M PBS pH 7.2,0.05% Tween<sup>®</sup> 20 four times
5. Add 100  $\mu$ l /well of TMB substrate solution and incubate 15 minutes at room temperature.
6. Stop the substrate reaction by adding 100  $\mu$ l /well of sulphuric acid 1 N and read the optical density values at 450 nm

The data show that a plateau has got starting with an IgG concentration of 8.0 $\mu$ g/ml.

This concentration means the well binding capacity we can express as :

- $\mu$ g/well = 0.800 ( 800 ng/well )
- pmol/well= 5.3 ( this result is calculated considering the IgG M.W. = 150.000 )

Figure 1



## TECHNICAL NOTES N. 27    sensitivity test

1. Add 100  $\mu$ l of different concentrations of human biotinylated IgG ( from 1.56 to 100 ng/ml ) to the wells of Protein G coated plate and incubate for 30 minutes at room temperature.
2. Empty the wells and wash with 0.1 M PBS pH 7.2,0.05% Tween<sup>®</sup> 20 four times.
3. Add 100  $\mu$ l /well of Streptavidin-HRP ( BioSpa product code SB01-61 , diluted 1:20.000 ) and incubate for 30 minutes at room temperature.
4. Empty the wells and wash with 0.1 M PBS pH 7.2,0.05% Tween<sup>®</sup> 20 four times.
5. Add 100  $\mu$ l /well of TMB substrate solution and incubate 15 minutes at room temperature.
6. Stop the substrate reaction by adding 100  $\mu$ l /well of sulphuric acid 1 N and read the optical density values at 450 nm

The microplate sensitivity was calculated as the lowest biotinylated IgG concentration higher than the mean optical density plus 5 S.D. of 0 ng/ml biotinylated IgG concentration.

Our experiment gave the following results :

- 0 ng/ml biotinylated IgG optical density mean ( coming from 8 replicates ) = 0.133
- standard deviation = 0.012
- mean + 5 S.D. = 0.193
- sensitivity = 0.056 ng/well of human IgG

sensitivity of protein G coated plate




Figure 2